


Deltalab 0.4ml Microtubes, Beckman Type

Product Code: 900007

Brand: Deltalab

Manufactured in polyethylene, Deltalab's 0.4ml Microtubes, Beckman Type can withstand centrifugation up to 11,000xg and can be used at temperatures from -10°C to 130°C.



Specifications

Dimensions (Ø x h, mm)	5.50 x 47.50
Description	Beckman® type microtube, 0.4ml
Quantity	1000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com