

Deltalab 0.7ml Cobas Bio HDPE Sample Cups

Product Code: 900008

Brand: Deltalab

Deltalab's 0.7ml Cobas Bio HDPE Sample Cups provide excellent chemical and temperature resistance.

Specifications

Volume (ml)	0.7
Type	Cobas Bio
External Diameter (mm)	7.65
Height (mm)	35.5
Material	HDPE
Quantity	1000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com