

Deltalab 2.5ml Hitachi Sample Cups

Product Code: 900020

Brand: Deltalab

Made of polystyrene, Deltalab's 2.5ml Hitachi Sample Cups are designed for use with Hitachi analysers.

Specifications

Volume (ml)	2.5
Type	Hitachi
External Diameter (mm)	16.75
Height (mm)	38
Material	Polystyrene
Sterility	Non-Sterile
Quantity	1000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com