

## Deltalab 0.5ml Gemsaec Sample Cups

**Product Code:** 900024

**Brand:** Deltalab

Deltalab's 0.5ml Gemsaec Sample Cups are made of polypropylene with polystyrene caps and are designed to fit popular analysers.

### Specifications

Volume (ml)	0.5
Type	Gemsaec, Kone Lab 20
External Diameter (mm)	13.55
Height (mm)	24.50
Material(s)	Polystyrene
Sterility	Non-Sterile
Quantity	1000



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)