

Deltalab Individually Wrapped Sterile Serological Pipettes (1ml) - Tape Ended

Product Code: 900030R

Brand: Deltalab


Made of clear polystyrene for perfect measurements, Deltalab's Individually Wrapped Sterile Serological Pipettes (1ml) - Tape Ended can be used for various laboratory applications (e.g. tissue culture). Each pipette is sterilised by radiation and made up of one piece for optimum precision.

Features:

- No liquid retention: as tip is not sealed, liquid is not held at the joint. Total dispensing guaranteed.
- Non-pyrogenic, non-cytotoxic (ISO 10993-5) and non-haemolytic (ISO 10993-4). DNase and RNase free.
- Include a coloured cotton plug for better identification.

Specifications:

Capacity (ml)	1
Cotton Colour	Yellow
Tip	Tape End
Graduation (ml)	0/0.9
Negative Graduation (ml)	Until -0.3
Total Capacity (ml)	1.3
Material	Polystyrene
Sterility	Sterile
Quantity	1000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com