

## Deltalab Bulk Wrapped Sterile Serological Pipettes (1ml) - Tape Ended

**Product Code:** 900031R

**Brand:** Deltalab


Made of clear polystyrene for perfect measurements, Deltalab's Bulk Wrapped Sterile Serological Pipettes (1ml) - Tape Ended can be used for various laboratory applications (e.g. tissue culture). Each pipette is sterilised by radiation and made up of one piece for optimum precision.

### Features:

- No liquid retention: as tip is not sealed, liquid is not held at the joint. Total dispensing guaranteed.
- Non-pyrogenic, non-cytotoxic (ISO 10993-5) and non-haemolytic (ISO 10993-4). DNase and RNase free.
- Include a coloured cotton plug for better identification.

### Specifications:

<b>Capacity (ml)</b>	1
<b>Cotton Colour</b>	Yellow
<b>Tip</b>	Tape End
<b>Graduation (ml)</b>	0/0.9
<b>Negative Graduation (ml)</b>	Until -0.3
<b>Total Capacity (ml)</b>	1.3
<b>Material</b>	Polystyrene
<b>Sterility</b>	Sterile
<b>Quantity</b>	1000



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)