

Deltalab Individually Wrapped Sterile Serological Pipettes (2ml)

Product Code: 900032R

Brand: Deltalab

Sterilised by radiation, Deltalab's Individually Wrapped Sterile Serological Pipettes (2ml) are made up of one piece for optimum precision.

Features

- No liquid retention: as tip is not sealed, liquid is not held at the joint. Total dispensing guaranteed
- Non-pyrogenic, non-cytotoxic (ISO 10993-5) and non-haemolytic (ISO 10993-4). DNase and RNase free
- Include a coloured cotton plug for better identification



Specifications

| | |
|---------------------------------|-------------|
| Capacity (ml) | 2 |
| Cotton Colour | Green |
| Tip | Tape End |
| Graduation (ml) | 0/1.8 |
| Negative Graduation (ml) | Until -0.4 |
| Total Capacity (ml) | 2.4 |
| Material | Polystyrene |
| Sterility | Sterile |
| Quantity | 800 |

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com