

Deltalab Bulk Wrapped Sterile Serological Pipettes (2ml)

Product Code: 900033R

Brand: Deltalab

Sterilised by radiation, Deltalab's Bulk Wrapped Sterile Serological Pipettes (2ml) are made up of one piece for optimum precision.

Features

- No liquid retention: as tip is not sealed, liquid is not held at the joint. Total dispensing guaranteed
- Non-pyrogenic, non-cytotoxic (ISO 10993-5) and non-haemolytic (ISO 10993-4). DNase and RNase free
- Include a coloured cotton plug for better identification.



Specifications

Capacity (ml)	2
Cotton Colour	Green
Tip	Tape End
Graduation (ml)	0/1.8
Negative Graduation (ml)	Until -0.4
Total Capacity (ml)	2.4
Material	Polystyrene
Sterility	Sterile
Quantity	1000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com