

Deltalab Individually Wrapped Tape Ended Sterile Serological Pipettes (10ml)

Product Code: 900036R

Brand: Deltalab

Execute fast, precise fluid dispensing for laboratory applications with Deltalab's Individually Wrapped Tape Ended Sterile Serological Pipettes (10ml).


Sterilised by radiation, each pipette is made up of one piece for optimum precision.

Features

- **No liquid retention:** as tip is not sealed, liquid is not held at the joint. Total dispensing guaranteed
- Non-pyrogenic, non-cytotoxic (ISO 10993-5) and non-haemolytic (ISO 10993-4). DNase and RNase free
- Include a coloured cotton plug for better identification

Specifications

Capacity (ml)	10
Cotton Colour	Orange
Tip	Tape End
Graduation (ml)	0.9
Negative Graduation (ml)	Until -3
Total Capacity (ml)	13
Material	Polystyrene
Sterility	Sterile
Quantity	200



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com