

## Deltalab Bulk Wrapped Sterile Serological Pipettes (5ml)

**Product Code:** 900038R

**Brand:** Deltalab

Execute fast, precise fluid dispensing for laboratory applications with Deltalab's Bulk Wrapped Sterile Serological Pipettes (5ml).

Sterilised by radiation, each pipette is made up of one piece for optimum precision.

### Features

- **No liquid retention:** as tip is not sealed, liquid is not held at the joint. Total dispensing guaranteed
- Non-pyrogenic, non-cytotoxic (ISO 10993-5) and non-haemolytic (ISO 10993-4). DNase and RNase free
- Include a coloured cotton plug for better identification



### Specifications

<b>Capacity (ml)</b>	5
<b>Cotton Colour</b>	Purple
<b>Tip</b>	Tape End
<b>Graduation (ml)</b>	0/4
<b>Negative Graduation (ml)</b>	Until -3
<b>Total Capacity (ml)</b>	8
<b>Material</b>	Polystyrene
<b>Sterility</b>	Sterile
<b>Quantity</b>	500

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)