

Deltalab Individually Wrapped Sterile Serological Pipettes (25ml)

Product Code: 900041R

Brand: Deltalab

Made of crystal polystyrene with tapered ends, Deltalab's Individually Wrapped Sterile Serological Pipettes (25ml) feature neat, easy-to-read graduations and can be used for various laboratory applications (e.g. tissue culture).

Features

- No liquid retention: as tip is not sealed, liquid is not held at the joint. Total dispensing guaranteed
- Non-pyrogenic, non-cytotoxic (ISO 10993-5) and non-haemolytic (ISO 10993-4). DNase and RNase free
- Include a coloured cotton plug for better identification



Specifications

Capacity (ml)	25
Cotton Colour	Dark Orange
Tip	Tape End
Graduation (ml)	0/23
Negative Graduation (ml)	Until -10
Total Capacity (ml)	35
Material	Polystyrene
Sterility	Sterile
Quantity	200

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com