


Dotalab Single Use Straw-Like Pipettes

Product Code: 900050

Brand: Dotalab

Suitable for automatic pipettes, Dotalab's Single Use Straw-Like Pipettes are made of polypropylene and are sterilised by radiation.



Specifications:

Capacity (ml)	1
Length (mm)	190
Case Quantity	40 x 25

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com