

Deltalab 20ml Scintillation Vials

Product Code: 900100

Brand: Deltalab

Deltalab's 20ml Scintillation Vials are single-use vials made of high density polyethylene designed to minimize solvent losses.

Their screw caps ensure a leakproof seal.

These vials are compatible with most liquid scintillation containers available on the market.



Specifications

Volume (ml)	20
Dimensions (Ø x h)	26.5 x 58.5mm
Quantity	1000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com