

## Deltalab 4ml Minivials for 20ml Scintillation Vials

**Product Code:** 900101

**Brand:** Deltalab


Deltalab's 4ml Minivials for 2ml Scintillation Vials are single-use vials made of high density polyethylene designed to minimize solvent losses.

Their screw caps ensure a leakproof seal.

These vials are compatible with most liquid scintillation containers available on the market.

### Specifications

<b>Volume (ml)</b>	4
<b>Dimensions (Ø x h)</b>	13.71 x 53.15mm
<b>Quantity</b>	2000



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)