

Deltalab 0.5ml Microtubes, Vitraton Type

Product Code: 900117

Brand: Deltalab

Made of autoclavable polypropylene, Deltalab's 0.5ml Microtubes, Vitraton Type are ideal for little sample volumes, microtechniques, microsedimentation, etc.

They are suitable for transport, storage, and freezing of biological samples.

They have a watertight closing system.


Their caps can be easily opened and closed with one hand. The rim in the inner part of the caps assure an optimum watertightness.

They are easily pierceable and will withstand centrifugation up to 11,000xg.

They can resist temperatures from -10°C to 130°C.

Specifications

Dimensions (Ø x h)	7.50 x 31.25mm (capped)
Volume (ml)	0.5
Quantity	1000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com