


Deltalab 0.5ml Microtubes, Vitraton Type (Graduated)

Product Code: 900118

Brand: Deltalab

Manufactured in autoclavable polypropylene, Deltalab's 0.5ml Microtubes, Vitraton Type (Graduated) will withstand centrifugation up to 14,000xg.

Their caps can be easily opened and closed with one hand. The rim in the inner part of the caps assure an optimum watertightness.



Specifications

Dimensions (Ø x h)	7.87 x 31.8mm (capped)
Volume (ml)	0.5
Quantity	1000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com