

## Deltalab 4ml Flat Bottom Screw Threaded Vials

**Product Code:** 900204

**Brand:** Deltalab

Made of ultra clear borosilicate glass, Deltalab's 4ml Flat Bottom Screw Threaded Vials are fully autoclavable and supplied in trays individually shrink-wrapped.



### Specifications

<b>Thickness (mm)</b>	1
<b>Dimensions (mm)</b>	15 x 45
<b>Volume (ml)</b>	4
<b>Screw</b>	13 - 425
<b>Quantity</b>	100

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)