

Deltalab 5ml Flat Bottom Screw Threaded Vials

Product Code: 900205

Brand: Deltalab

Made of ultra clear borosilicate glass, Deltalab's 5ml Flat Bottom Screw Threaded Vials are fully autoclavable and supplied in trays individually shrink-wrapped.



Specifications

Thickness (mm)	1
Dimensions (mm)	18 x 50
Volume (ml)	5
Screw	15 - 425
Quantity	250

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com