

Deltalab 20ml Flat Bottom Screw Threaded Vials

Product Code: 900220

Brand: Deltalab

Made of ultra clear borosilicate glass, Deltalab's 20ml Flat Bottom Screw Threaded Vials are fully autoclavable and supplied in trays individually shrink-wrapped.

Specifications

Thickness (mm)	1
Dimensions (mm)	28 x 65
Volume (ml)	20
Screw	22 - 400
Quantity	300



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com