


Deltalab 4ml Round Bottom Borosilicate Glass Tubes

Product Code: 901075

Brand: Deltalab

The high quality of Deltalab's 4ml Round Bottom Borosilicate Glass Tubes are reflected in the uniformity of their wall thickness and of their diameter and height dimensions.

They are supplied in small quantities per case for a more convenient use in laboratories.



Specifications

Total Capacity (ml)	4
Internal Diameter, tube (mm)	8.20
External Diameter, tube (mm)	9.75
Height (mm)	75
Thickness (mm)	0.60
Quantity	1000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com