

Deltalab Single Use Polystyrene Racks (13mm)

Product Code: 901410

Brand: Deltalab

Deltalab's Single Use Polystyrene Racks (13mm) can be used for storing large quantities of test tubes. Each laboratory rack features deep holes to provide support for the tubes.



Specifications

Dimensions (mm)	160 x 170 x 25
Hole Diameter (mm)	13
Number of Holes	10 x 10
Material	Polystyrene
Sterility	Non-Sterile
Quantity	25

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com