

Ceti 1.3 Mpx. Digital USB Camera

Product Code: 9554.1200M

Brand: Ceti

The Ceti 1.3 Mpx. Digital USB Camera comes with eye piece adaptors to fit into an eye piece tube once it has been removed.

This microscope camera's adaptors can be used with different diameter eye pieces, which allows simple software to effectively capture images.



Specifications:

Sensor Size	1/3"
Pixel Size	8µm x 8µm
Sensitivity (v/lux - sec.)	1.8
Max. Resolution (hardware)	1280 x 1024 1.3 Mpx.
Exposure	ERS (Electronic Rolling Snap)
White Balance	Automatic/Advanced Manual
Auto-Exposure	Automatic/Advanced Manual
Image Output	USB 2.0
Programmable Control	Image size, exposure time, brightness, gain, gamma correction, colour, contrast, saturation
Accessories Included	Two adaptors (23.2 to 30, 23.2 to 30.50, one each) TS-M1 stage micrometer available at extra cost
Output Data	RGB24, Direct Show, and Twain Interface
Software	Driver: support Direct Show and Twain Image software: MiniSee

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com