

Ceti Digi-Pad 5 Mpx. Microscope Tablet Camera (With 9.7" Screen)

Product Code: 9556.5600M

Brand: Ceti

Ceti's Digi-Pad 5 Mpx. Microscope Tablet Camera (With 9.7" Screen) is specially designed for microscopy applications. It comes with a 9.7" colour LCD high-definition screen and is equipped with special software that supports, previews, photographs, video measurement capabilities, annotations, and saving.

This modern microscope camera also has a C-mount connector, which can be connected to any CETI microscope with an adaptor.

Microscope not included.

Features

- 9.7 inch colour LCD high-definition screen with true colour reproduction
- With built-in special microscopy software, the camera runs smoothly to give a perfect user experience
- Crystal display, RK3066 dual-core 1.5Hz CPU
- Touch-screen operation, can be used as a small computer. Keyboard and mouse can also be used to operate
- Includes Wi-Fi and Bluetooth wireless networking for easy file exchange



Specifications

CPU	RK3066 dual-core 1.5Hz
Display	9.7 inch (1024 x 768) colour LCD touch-screen
Storage	External SD memory card - up to 32GB
	ROM 2GB
	RAM 1GB DDR3
Camera	15fps at 1024 x 768
	1/2.5 (4:3) 5 mpx. CMOS
Network	Wi-Fi
	Bluetooth
I/O	HDMI Output
	On/off Button
	DC5V Power
	SD Memory Slot
	USB/Mini-USB (OTG) supports USB keyboard and mouse
Mechanical	LCD screen 360° horizontal/90° vertical
Operating Software	Android 4.2 RAM 2GB
Software	Built-in microscopy software

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com