

## Ceti 11.6" HD Monitor for WFHD5C Camera

**Product Code:** 9556.5800M

**Brand:** Ceti

Ceti's 11.6 HD Monitor for WFHD5C Camera features a Panasonic IPS LCD (Super TFT) to guarantee a wide viewing angle and high resolution. The LCD panel has true 1080P resolution, a high contrast ratio (up to 1000:1), LED backlight with 50,000 hours long life time and an 11.6" active area.

The screen is protected by a high impact ABS bezel and the back of the monitor includes a screw fitting, allowing the device to be mounted on the camera. A twist and tilt fitting allows the user to adjust the viewing angle.



# Specifications

<b>LCD Panel</b>		Panasonic IPS LCD screen (Super TFT)
<b>Input Video Format</b>		HDMI
<b>Native Resolution</b>		1920 x 1080
<b>Display Type</b>		16:9 Ratio, 11.6" Active Matrix Super TFT LCD
<b>Typical Contrast Ratio</b>		1000:1
<b>Colours</b>		16.7 million
<b>Viewing Angle (L/R/U/D)</b>		IPS Full View
<b>Active Display Area (W x H)</b>		258 x 145mm
<b>Pixel Pitch (W x H)</b>		0.134 x 0.134mm
<b>Brightness</b>		350cd/sq.m; 400cd/sq.m (optional)
<b>Backlight</b>		LED backlight: 50,000 hours
<b>Housing</b>	<b>Colour</b>	Black
	<b>Dimensions (L x H x W)</b>	281 x 179 x 15.6mm
	<b>Weight (g)</b>	400
<b>Operating Environment</b>	<b>Operating Temperature (°C)</b>	-15 to +55
	<b>Operating Humidity (%RH)</b>	10 to 90
	<b>Storage Humidity (%RH)</b>	5 to 90
	<b>Synchronization Range</b>	30 - 80 KHz Horizontal, 55 - 75 KHz Vertical
	<b>Power Supply</b>	AC 110V - 220V/DC12V (1A)
	<b>Power Consumption</b>	Max. 12W

Unit 3, Tower Business Park  
 Warpsgrove Lane  
 Chalgrove  
 Oxfordshire  
 OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)