


Deltalab Stainless Steel Rack (for 50 x 18mm Tubes)

Product Code: A-210

Brand: Deltalab

Manufactured from high-quality stainless steel, Deltalab's Stainless Steel Rack (for 50 x 18mm Tubes) can hold up to 50 tubes with a diameter of 18mm. It has holes with a diameter of 19mm that can firmly keep the tubes in the correct position.



Specifications:

Dimensions (l x w x h)	210 x 110 x 60mm
Containers Capacity	50 (10 x 5)
Material	Stainless Steel
Sterility	Non-Sterile
Quantity	1

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com