

## Deltalab Stainless Steel Rack (for 50 x 18mm Tubes)

**Product Code:** A-210

**Brand:** Deltalab

Manufactured from high-quality stainless steel, Deltalab's Stainless Steel Rack (for 50 x 18mm Tubes) can hold up to 50 tubes with a diameter of 18mm. It has holes with a diameter of 19mm that can firmly keep the tubes in the correct position.

### Specifications:

<b>Dimensions (l x w x h)</b>	210 x 110 x 60mm
<b>Containers Capacity</b>	50 (10 x 5)
<b>Material</b>	Stainless Steel
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	1

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)