

Bigneat Chemcap Clearview Ductless Fume Cabinet (BC I206)

Product Code: BIGBCI206

Brand: Bigneat

Bigneat's Chemcap Clearview Ductless Fume Cabinet (BC I206) operates with excellent airflow control, highly efficient carbon filtration, is energy efficient and easy to use. View of the working area is an unobstructed all-round 360° and this makes the cabinet comfortable to use for long periods.

Supplied complete with 2x pre-filters, 2x carbon filters of pre-determined type and 2x HEPA filters on request (double filter stack and double fan system).

Features:

- Keypad control panel with touch button controls and service due indicator.
- Audible electronic low airflow alarm.
- Clearview all-round vision.
- Easy clean - all materials chemically resistant.
- Energy efficient fan system.
- Energy efficient LED lighting.
- Lockable front access window.
- Pre-filter can be bagged under negative pressure.
- Quiet operation.



Specifications:

Dimensions (w x d x h)	Internal	1170 x 520 x 810mm
	External	1200 x 570 x 1145mm
Access Style		Trapezium
Face Velocity		0.5 m/s
Fan System Rated		104W
Average Power Consumption		80W
Sound Level (dBA)		< 55

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com