

Medline BECELEC 2 Electric Burner (2 Pin EU)

Product Code: BEC2EXSET

Brand: Medline Scientific

The Medline BECELEC 2 Electric Burner (2 PIN EU) offer laboratories a unique safe replacement option to conventional gas Bunsen burners. The high power heater combined with a reflective internal stainless steel bowl produces a concentrated stream of heat comparable with that of existing gas-powered Bunsen burners.

Supplied with reflective internal stainless steel bowl, wire safety guard, retort stand, support rod, boss, and clamps. One litre gel flame optional.



Specifications:

Temperature Range (°C)	37 to 720
Diameter of Heating Area (mm)	Ø 130
Power	500 Watt
Safety Indicator	LED
Power Supply	220V AC
Plug Type	2 pin EU
Dimensions (w x d x h)	31 x 20 x 11cm
Weight (kg)	1.7
Accessories	Stainless steel plate, flame gel, support rod, boss head, and clamp

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com