

Neptune Racked S3 Polymer Filter Tips (1000 - 1250µl)

Product Code: BT1250

Brand: Neptune Scientific

Made of virgin polypropylene, Neptune's Racked S3 Polymer Filter Tips (1000 - 1250µl) are pre-sterilised and S3 polymer tips are used to stop samples being held up and carry out a precise sample draw and dispense.



Specifications:

Graduation Marks (µl)	1000 - 1250
Tip Composition	Virgin Polypropylene
Tip Type	Natural Polypropylene and S ³
Filter Material	High Density Polyethylene Filter
Offered in Pre-Sterile Format?	Yes
Configuration	Racked
Autoclavable?	No
Quantity	8 x 96

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com