

## Neptune Racked S3 Polymer Filter Tips (200µl)

**Product Code:** BT200

**Brand:** Neptune Scientific

Made of virgin polypropylene, Neptune's Racked S3 Polymer Filter Tips (200µl) are pre-sterilised and S3 polymer tips are used to stop samples being held up and carry out a precise sample draw and dispense.



### Specifications:

<b>Graduation Marks (µl)</b>	200
<b>Tip Composition</b>	Virgin Polypropylene
<b>Tip Type</b>	S <sup>3</sup>
<b>Filter Material</b>	High Density Polyethylene Filter
<b>Offered in Pre-Sterile Format?</b>	Yes
<b>Configuration</b>	Racked
<b>Autoclavable?</b>	No
<b>Quantity</b>	10 x 96

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)