

Jeio Tech Dual Water Bath (BW-0505H)

Product Code: AAH47115K

Brand: Jeio Tech

The Jeio Tech Dual Water Baths (BW-0505H) is an economical, digital heating water bath that combine two baths into one.

This high-quality, robust laboratory water bath is used for the direct immersion of samples and provides great temperature stability for various laboratory applications.

Features:

Performance:

- Temperature range from ambient +5°C to 100°C.
- Microprocessor PID control/automatic tuning/calibration.
- Powerful circulation system using our patented agitator.

Convenience:

- Two independently controlled reservoirs merged into a single unit for space saving.
- Corrosion-resistant stainless steel construction allowing other bath fluids as well as water and silicone fluids.
- Digital LED displays with quick adjustable dial knob (0.1°C resolution).
- Low-profile design with sloped control panel for easier access and better safety. Also side handles for easy transport.

Safety:

- Isolated design of heater and sensor by a baffle plate to protect users and samples against direct contact.
- Complete safety protection system with warning alarms. Over-temperature, over-current protection.
- Keypad lock function prevents accidental changes.



Specifications:

Bath Volume (litres)	3.5 & 3.5	
Working Temperature Range (°C)	Amb. +5 to 100	
Temperature Stability (°C)	± 0.1	
Dimensions	Bath Opening/Depth (w x l, d)	240 x 136, 150mm per chamber
	External (w x d x h)	590 x 200 x 265mm
	Net Weight (kg)	9.7
Electrical Requirements (230V, 50/60Hz)	3 & 3A	

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com