

Jeio Tech 3.5 Litre Water Bath (BW-05B)

Product Code: AAH46012K

Brand: Jeio Tech

Constructed with a corrosion-resistant stainless steel tank, the Jeio Tech 3.5 Litre Water Bath (BW-05B) is an economical water bath that provides stable temperature control for samples.

This unit comes with an easy and fast sloped control panel, a safety alarm setting, and splash-proof keypad.

Features:

Performance

- Temperature range from ambient +7°C to 100°C.
- On/off temperature control.

Convenience

- Simple and quick control panel.
- Corrosion resistant stainless steel construction allowing other bath fluids as well as water and silicone fluids.
- Low-profile design with sloped control panel for easier access and better safety.
- Variety of accessories are available to meet application demands.

Safety

- Isolated design of heater and sensor by a baffle plate to protect users and samples against direct contact.
- Complete safety protection system with warning alarms. Over-temperature, over-current protection.
- Splash-proof keypad.



Specifications:

Bath Volume (litres)		3.5
Working Temperature Range (°C)		Amb. +7 to 100
Temperature Stability (°C)		± 0.5
Dimensions	Bath Opening/Depth (w x l, d)	240 x 136, 115mm
	External (w x d x h)	303 x 242 x 245mm
	Net Weight (kg)	6.5
Electrical Requirements (230V)		50Hz, 3A

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com