

Jeio Tech Dual Water Bath (BW-I020H)

Product Code: AAH47145K

Brand: Jeio Tech

The Jeio Tech Dual Water Bath (BW-I020H) is an economical, digital heating water bath that combine two baths into one.

This high-quality, robust laboratory water bath is used for the direct immersion of samples and provides great temperature stability for various laboratory applications.

Features:

Performance:

- Temperature range from ambient +5°C to 100°C.
- Microprocessor PID control/automatic tuning/calibration.
- Powerful circulation system using our patented agitator.

Convenience:

- Two independently controlled reservoirs merged into a single unit for space saving.
- Corrosion-resistant stainless steel construction allowing other bath fluids as well as water and silicone fluids.
- Digital LED displays with quick adjustable dial knob (0.1°C resolution).
- Low-profile design with sloped control panel for easier access and better safety. Also side handles for easy transport.

Safety:

- Isolated design of heater and sensor by a baffle plate to protect users and samples against direct contact.
- Complete safety protection system with warning alarms. Over-temperature, over-current protection.
- Keypad lock function prevents accidental changes.



Specifications:

Bath Volume (litres)		11.5 & 20
Working Temperature Range (°C)		Amb. +5 to 100
Temperature Stability (°C)		± 0.1
Dimensions	Bath Opening/Depth (w x l, d)	300 x 240, 200mm and 498 x 300, 200mm
	External (w x d x h)	844 x 372 x 318mm
	Net Weight (kg)	17.6
Electrical Requirements (230V, 50/60Hz)		4.3A & 8.7A

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com