

Jeio Tech 11.5 Litre Stirred Water Bath (BW-10H)

Product Code: AAH45125K

Brand: Jeio Tech

The Jeio Tech 11.5 Litre Stirred Water Bath (BW-10H) features a microprocessor PID control with automatic tuning and calibration plus digital LED displays with a quick adjustable dial knob. The powerful circulation system, using the patented agitator, ensures temperature accuracy of $\pm 0.1^{\circ}\text{C}$.

Corrosion resistant stainless steel construction allows the use of other bath fluids as well as water and silicone oil. Side handles are included for easy transport. Safety is ensured with over-temperature and over-current protection. A keypad lock prevents accidental changes.

Features

Performance

- Temperature range from ambient $+5^{\circ}\text{C}$ to 100°C
- Microprocessor PID control/automatic tuning/calibration
- Powerful circulation system using our patented agitator

Convenience

- Corrosion resistant stainless steel construction allowing other bath fluids as well as water and silicone fluids
- Digital LED displays with quick adjustable dial knob (0.1°C resolution)
- Low-profile design with sloped control panel for easier access and better safety. Also side handles for easy transport

Safety

- Isolated design of heater and sensor by a baffle plate to protect users and samples against direct contact
- Complete safety protection system with warning alarms. Over-temperature, over-current protection
- Keypad lock function prevents accidental changes



Specifications

Bath Volume (L)	11.5
Temperature Range (°C)	Amb. +5 to 100
Temperature Accuracy (°C)	± 0.1
Bath Opening/Depth (W x L, D)	300 x 240, 200mm
Dimensions (W x D x H)	364 x 316 x 318mm
Weight (kg)	9.5

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com