


Jeio Tech 5 to 100rpm Rocking Shaker (CRS-350)

Product Code: AAH391115K

Brand: Jeio Tech

The Jeio Tech 5 to 100rpm Rocking Shaker (CRS-350) comes with an adjustable tilt angle and precise digital speed control, providing a desirable rocking motion.

Three different detachable trays for this laboratory shaker can be used to efficiently transfer well-plates and test tubes in various sizes.



Features

Performance

- Microprocessor PID control
- Digital speed control providing gentle to vigorous motion (5 to 100rpm)
- Easy to adjust tilt angle (7 to 13°)
- Smooth, quiet, uniform, and yet powerful shaking motion

Convenience

- Vivid and colourful fluorescent display (VFD) with touch-sensitive keypad
- One-touch timer for time-dependent shaking operation (1 min. to 99 hr. 59 min.)
- Three types of removable tray for easy transport of well-plates, test tubes, etc. (standard dimpled mat, and optional rubber mat and universal harness)
- Variety of accessories are available to meet application demands. Accommodates to be used in a variety of test tube size, beakers, vials, micro well-plates, culture dishes, plastic/glass trays

Safety

- Self-diagnostic function identifying errors
- Its polypropylene main body is highly chemical-resistant and easy-to-clean
- Over-current protection

Specifications

Control System	PID feedback control
Display	VFD
Motion	Rocking
Speed Range (rpm)	5 to 100
Tilt Angle	7 to 13°
Max. Load (kg)	5
Tray Dimensions (W x D)	276 x 276mm
Dimensions (W x D x H)	249 x 325 x 138mm
Weight (kg)	3.7

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com