

## Deltalab EUROTUBO Cut 90° Slides (No Matt Area)

**Product Code:** D100001

**Brand:** Deltalab

Ideal for routine analysis, Deltalab's EUROTUBO Cut 90° Slides (No Matt Area) are made of clear pre-washed glass. Displayed in vacuum flow-packed cases of 50 units, these microscope slides are protected from dirtiness and include a bag of silica salt.



### Specifications:

<b>Edges</b>	Cut 90°
<b>Corners</b>	Non-clipped
<b>Matt Area</b>	No
<b>Standard Dimensions</b>	26 x 76mm
<b>Thickness</b>	1.1mm
<b>Expiry Date</b>	24 months after manufacturing date
<b>Material</b>	Glass
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	50

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)