


Deltalab EUROTUBO Grounded 45° Slides (No Matt Area)

Product Code: D100002

Brand: Deltalab

Ideal for routine analysis, Deltalab's EUROTUBO Grounded 45° Slides (No Matt Area) are made of clear pre-washed glass. Displayed in vacuum flow-packed cases of 50 units, these microscope slides are protected from dirtiness and include a bag of silica salt.



Specifications:

Edges	Grounded 45°
Corners	Clipped
Matt Area	No
Standard Dimensions	26 x 76mm
Thickness	1.1mm
Expiry Date	24 months after manufacturing date
Material	Glass
Sterility	Non-Sterile
Quantity	50

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com