

## Deltalab 0.5ml Skirted Sterile Screw Cap Microtubes with Caps (Non-Graduated)

**Product Code:** 409115/4

**Brand:** Deltalab

Made of medical-grade and transparent polypropylene, Deltalab's 0.5ml Skirted Sterile Screw Cap Microtubes with Caps (Non-Graduated) are supplied in bags of 50 units.

They are suitable for autoclave, liquid nitrogen (gas) and boiling processes and are perfect for long-term storage.


They can withstand temperatures from -190°C and are DNase, RNase, DNA, and PCR inhibitors free.

The caps embody a non-reactive propylene o-ring.

These microtubes incorporate an external grip for easy handling.

### Specifications

<b>Description</b>	0.5ml skirted
<b>Graduation</b>	No
<b>Material</b>	Polypropylene
<b>Sterility</b>	Sterile
<b>Quantity</b>	2500



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)