

## Deltalab 1.5ml Non-Skirted Sterile Screw Cap Microtubes with Caps (Graduated)

**Product Code:** 409115/2

**Brand:** Deltalab

Made of medical-grade and transparent polypropylene, Deltalab's 1.5ml Non-Skirted Sterile Screw Cap Microtubes with Caps (Graduated) are supplied in bags of 50 units.

They are suitable for autoclave, liquid nitrogen (gas) and boiling processes and are perfect for long-term storage.


They can withstand temperatures from -190°C and are DNase, RNase, DNA, and PCR inhibitors free.

The caps embody a non-reactive polypropylene o-ring.

These microtubes feature a frosted area for writing.

### Specifications

<b>Description</b>	1.5ml non-skirted
<b>Graduation</b>	Yes
<b>Material</b>	Polypropylene
<b>Sterility</b>	Sterile
<b>Quantity</b>	2500



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)