

Deltalab 1.5ml Transparent Polypropylene Screw Thread Microtubes (Skirted)

Product Code: 409111/5

Brand: Deltalab


Deltalab's 1.5ml Transparent Polypropylene Screw Thread Microtubes (Skirted) are suitable for use in liquid nitrogen, autoclave, and for boiling applications. They can be used at temperatures down to -190°C .

They are DNase, RNase, and PCR inhibitors free.

These microtubes are ideal for long-term storage.

They can withstand centrifugation up to 20,000xg.

Caps are supplied separately.



Specifications

Description	1.5ml (skirted)
Dimensions (mm)	10.3 x 44.5
Graduated?	No
Quantity	500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com