

Deltalab 2ml Non-Skirted Sterile Screw Cap Microtubes with Caps (Graduated)

Product Code: 409115/3

Brand: Deltalab

Made of medical-grade and transparent polypropylene, Deltalab's 2ml Non-Skirted Sterile Screw Cap Microtubes with Caps (Graduated) are supplied in bags of 50 units.

They are suitable for autoclave, liquid nitrogen (gas) and boiling processes and are perfect for long-term storage.

They can withstand temperatures from -190°C and are DNase, RNase, DNA, and PCR inhibitors free.

The caps embody a non-reactive propylene o-ring.

These microtubes feature a frosted area for writing.

Specifications

Description	2ml non-skirted
Graduation	Yes
Material	Polypropylene
Sterility	Sterile
Quantity	2500



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com

Features: