

## Deltalab 2ml Transparent Polypropylene Screw Thread Microtubes (Not Skirted)

**Product Code:** 409111/3

**Brand:** Deltalab


Deltalab's 2ml Transparent Polypropylene Screw Thread Microtubes (Not Skirted) are suitable for use in liquid nitrogen, autoclave, and for boiling applications. They can be used at temperatures down to -190°C.

They are DNase, RNase, and PCR inhibitors free.

These microtubes are ideal for long-term storage.

They can withstand centrifugation up to 20,000xg.

Caps are supplied separately.



### Specifications

<b>Description</b>	2ml (not skirted)
<b>Dimensions (mm)</b>	10.3 x 44.5
<b>Graduated?</b>	Yes
<b>Quantity</b>	500

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)