

## Deltalab 2ml Transparent Polypropylene Screw Thread Microtubes (Skirted & Graduated)

**Product Code:** 409111/7

**Brand:** Deltalab


Deltalab's 2ml Transparent Polypropylene Screw Thread Microtubes (Skirted & Graduated) are suitable for use in liquid nitrogen, autoclave, and for boiling applications. They can be used at temperatures down to  $-190^{\circ}\text{C}$ .

They are DNase, RNase, and PCR inhibitors free.

These microtubes are ideal for long-term storage.

They can withstand centrifugation up to 20,000xg.

Caps are supplied separately.



### Specifications

|                        |               |
|------------------------|---------------|
| <b>Description</b>     | 2ml (skirted) |
| <b>Dimensions (mm)</b> | 10.3 x 44.5   |
| <b>Graduated?</b>      | Yes           |
| <b>Quantity</b>        | 500           |

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)