


Deltalab Blue Cryovial Caps

Product Code: 409110/1

Brand: Deltalab

Made of polypropylene, Deltalab's Blue Cryovial Caps feature a long skirt and special thread design allowing them to be removed or sealed with a single turn and assuring a better protection of the content. A silicone washer ensures a positive leakproof seal at all temperatures.



Specifications:

Colour	Blue
Material	Polypropylene
Quantity	1000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com