

## Deltalab Blue STAR FROST Slides

**Product Code:** A100013

**Brand:** Deltalab


Made of high-definition pre-washed glass, and hydrophilic (epoxy), Deltalab's Blue STAR FROST Slides feature a colour-coded, 20mm marking area.

This special colour-printed area is resistant to all common chemicals used in laboratories, and may be marked with pencil, ballpoint pen, marker, and different laboratory printers (e.g. Leica XPS).

Each slide has ground edges at 45° (angle) and clipped corners.

### Specifications:

<b>Strip Colour(s)</b>	Blue
<b>Standard Size (mm)</b>	26 x 76
<b>Thickness (mm)</b>	1
<b>Expiry Date</b>	14 months after manufacturing date
<b>Material</b>	Glass
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	1000



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)