

Deltalab EUROTUBO Leak Proof Screw Cap Containers (1000ml)

Product Code: 202847

Brand: Deltalab

Specially designed for histology, Deltalab's EUROTUBO Leak Proof Screw Cap Containers (1000ml) provide double internal security closure.

Manufactured in translucent polypropylene with polyethylene ribbed yellow caps (supplied uncapped), they have moulded graduations in increments of 50ml and a ring in the upper part of the containers that prevent dripping when decanting liquids.

Their leak proof closure makes these rigid plastic containers ideal for liquid (including alcohols) or solid sample storage.



Specifications:

Capacity (ml)	1000
External Mouth (Ø mm)	111
Height (mm)	139
Material(s)	Polypropylene (container) Polyethylene (cap)
Sterility	Non-Sterile
Quantity	100

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com