

## Deltalab EUROTUBO Leak Proof Screw Cap Containers (500ml)

**Product Code:** 202846

**Brand:** Deltalab

Specially designed for histology, Deltalab's EUROTUBO Leak Proof Screw Cap Containers (500ml) provide double internal security closure. Manufactured in translucent polypropylene with polyethylene ribbed yellow caps (supplied uncapped), they have moulded graduations in increments of 50ml and a ring in the upper part of the containers that prevent dripping when decanting liquids.

Their leak proof closure makes these rigid plastic containers ideal for liquid (including alcohols) or solid sample storage.



### Specifications:

Capacity (ml)	500
External Mouth (Ø mm)	90
Height (mm)	110
Material(s)	Polypropylene (container) Polyethylene (cap)
Sterility	Non-Sterile
Quantity	200

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)