


## Deltalab Extra Long Pasteur Pipettes (23ml)

**Product Code:** 200060

**Brand:** Deltalab

Made of transparent polyethylene, Deltalab's Extra Long Pasteur Pipettes (23ml) carry out secure and simple sampling. They are suitable for tissue culture applications and each pipette features a built-in pipette bulb that helps prevent cross-contamination.

Each pipette is graduated each ml up to 5ml.



### Specifications:

|                                |              |
|--------------------------------|--------------|
| <b>Volumeterical Tolerance</b> | ± 10%        |
| <b>Length (mm)</b>             | 300          |
| <b>Stem (Ø mm)</b>             | 9            |
| <b>Total Capacity (ml)</b>     | 23           |
| <b>Bulb Draw (ml)</b>          | 7.1          |
| <b>Material</b>                | Polyethylene |
| <b>Sterility</b>               | Non-Sterile  |
| <b>Quantity</b>                | 100          |

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)