


Deltalab Extra Long Pasteur Pipettes (6ml)

Product Code: 210006

Brand: Deltalab

Made of transparent polyethylene, Deltalab's Extra Long Pasteur Pipettes (6ml) carry out secure and simple sampling. They are suitable for tissue culture applications and each pipette features a built-in pipette bulb that helps prevent cross-contamination.



Specifications:

Volumeterical Tolerance	± 10%
Length (mm)	225
Stem (Ø mm)	5
Total Capacity (ml)	6
Bulb Draw (ml)	2.3
Material	Polyethylene
Sterility	Non-Sterile
Quantity	400

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com