


Deltalab Filter Paper (50 x 50cm)

Product Code: I20010

Brand: Deltalab

Deltalab's Filter Paper (50 x 50cm) is made from high quality cellulose fibres to ensure a good resistance to humidity and a high absorption capacity. It can be used for the following applications: laboratory worktop protection, simple filtering operations for varied products, labware sterilisation, and pharmaceutical product preparation.



Specifications:

Sheet Size (cm)	50 x 50
Grammage	73g/m ²
Thickness (mm)	0.170 ± 0.02
Klemm Absorp. (mm)	89
Pore Diameter	± 30 - 40 microns
Material	Paper
Sterility	Non-Sterile
Quantity	100

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com